TendonCARE™
Ultrasound
 Therapy Massager
 Bio-838

CE Quality Certified
 Made in Hong Kong

Musculoskeletal Solution

- * **Speed up healing**
 Musculoskeletal Injuries
 Neuromuscular Disorders
 Bone Fractures
- * **Improves**
 Joints immobility
 Hamstring flexibility
 Blood Flow of Capillary
- * **Reduces**
 Resolution of Edema
 Softening Scar Tissue
- * **Power Penetration**
 1MHz 5cm↑ deep treatment
 2 W/cm² Therapy Power
- * **Functional Flexible**
 1~9 Mode digital display
 Presetable treatment time
 Continuous or Pulse mode
- * **Safe Treatment**
 Intelligent Monitoring
 LED and Buzzer Alert
 Power Diagnostic

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Those terms of disease do not use for advertising, please consult doctor)

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User manual and specification packed inside

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www.TendonCare.com

TendonCare Biophysics Sports Injuries Rehabilitation Products

TendonCARE™ (Bio-838) Ultrasound Therapy Massager

Mode Intensity

Applications:

Therapeutic ultrasound is one of the most common treatments used in the management of soft tissue lesions, which constitute the majority of rheumatic complaints.

Instructions of use: (In lack of ultrasound gel, ultrasound treatment effectiveness will be greatly reduced.)

1. Please consult your physiotherapist or doctor for correct setting on using this device.
2. Remove the treatment horn, if you feel burn or painful.
3. Connect the plug with your ultrasound device, before plugging on the power adaptor.
4. Preset the treatment mode and then time (see table on right)
5. Apply ultrasound gel on pained portion completely, and then gently keep it slowly swinging the horn around lesions
6. Do not keep the horn stationary on portion more than 30s. More than 60mins, sound alert will remind.
7. Ultrasound is a kind of sound, hard to feel vibration but a little warm. So, pay attention to patients, who do not feel enough feeling feedback. Horn is designed for working below 42°C. Device stops to work and LED indicator lights up if it is overrun.
8. Remove the power adaptor when the device is not in use. For every overrun, LED lights up, please turn off power at least 8min before restarting it.

Restrictions on the Use of Ultrasound Massager

- Over organs (stomach, liver, lungs, spleen, bowels, heart, kidneys, thyroid)
- Over the pelvis, abdomen and lumbar (lower back) region in pregnant or potentially pregnant females.
- On patients with pace-makers on patients with phlebitis, deep venous thromboses or patients with hemophilia
- Over mucous membranes (mouth, nose, rectum/anus, vagina)
- Over highly ultrasound-sensitive organs (eyes, ears, ovaries, testicles, brain, spinal cord)
- On patients with spinal bifida / bone with active infection / a bone fracture
- Over areas suspected to be cancerous / pre-cancerous
- Over de-sensitized (numb, hypoesthesia) areas of the skin (eg. diabetic neuropathy)
- Do not apply ultrasound gel on injured skin with wound
- Never use Ultrasound unit on children or under the age of 16 without the recommendation/supervision of a physician or therapist

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Warranty : 1year

Out of Warranty conditions :

- Water gets into device cavity
- Horn hits with hard materials
- Casing is opened unofficially
- Wrong electrical input

Specification

Frequency : 1MHz	Device Power : 12 W	Package : PP plastic box	Treatment operational Protection :
Pulse width : 2.5msec	Device Voltage input : 12V	Pack with: Ultrasound gel (450ml)	Horn temp. service : 42°C
Ultrasound Power : 0.03~2 W/cm ² (SATA)	Adaptor voltage : 220Vac / 110Vac	1pcs AC adaptor (220V / 110V)	Vibration detection Buzzer alert
Beam uniformity Ratio (BNR) : 4:1	Water resistance : No	User manual (this sheet)	Device overrun : LED alert
Effect Radiating Area (ERA) : 3	(can be cleaned by wet cloth)	Ultrasound device	Power connection : LED indicator
		Total weight : 1.5kg	Ultrasound output : LED indicator

Warnings: If any session of ultrasound therapy causes physical discomfort, stop using immediately. Prolonged exposure to maladjusted frequencies can cause serious tissue and neurological damage.

Attention: Correct setting of Ultrasound Therapy should not be painful and low temperature raised. Over setting would cause burning due to cavitation

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TendonCare Biophysics Co.Ltd